LABCONCO CORPORATION

8811 Prospect Avenue

Kansas City, MO 64132

(800) 821-5525, (816) 333-8811

(816) 363-0130 fax

labconco@labconco.com

Quick Start Installation Guide

**CentriVap® Centrifugal Concentrator Systems
Quick Guide to Installation**

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CentriVap® Centrifugal Concentrator Systems


Quick Guide to Installation

Note: This Quick Guide to Installation may not provide all the information for your installation. Refer to the User's Manual for details and precautions. Do not discard packaging until you have identified all parts that were shipped and the products operate properly.

Preparations

Step 1. Proper Electrical Supply


- 230V, 50 Hz models
 - 15-Amp dedicated single phase circuit. Additional appliances on the same circuit may cause the current to exceed the rating of the circuit breaker or fuse.
 - Frequency must agree with the serial tag rating. Improper frequency will damage the compressor.


UK plug

Step 2. Vacuum Pump Requirements

- For aqueous samples
 - 90 liters per minute free flow rate (5,2 m³ per hour)
 - 0,002 mBar ultimate vacuum
- For volatile solvent samples
 - 42 liters per minute (8,6 m³ per hour)
 - 200 mBar ultimate vacuum is adequate
 - PTFE
- Plug Requirements on vacuum pump power cord
 - 230V models - reverse IEC plug
- Hose barb connection to fit 13mm (0.5 inch) inside diameter vacuum hose


EU plug

Reverse IEC plug

Step 3. Vacuum Trap Requirements

- When using a mechanical pump that is not corrosion resistant, the CentriVap Cold Trap or equivalent trap must be used in the vacuum line to prevent damage from solvent vapors. When evaporating acids, it is advisable to include a soda lime acid secondary trap in the vacuum line. This adds additional protection for the pump



WARNING: It is recommended that the vacuum pump be located inside a fume hood or other laboratory ventilation device if hazardous solvents are used in the CentriVac. If this is not possible, the vacuum pump should have a hose attached to the exhaust port and the other end of the hose should be positioned inside the fume hood or ventilation device.


Installation

Step 4. Vacuum Pump Set Up (Rotary Vane pump shown)

- Remove protective caps from inlet and outlet ports.
- Some vacuum pumps have inlet fittings that screw into the inlet port. Make sure the O-ring is on the sealing surface. Insert fitting into the pump housing and tighten securely
- Some vacuum pumps have inlet fittings that clamp to the pump. Place the centering ring with filter screen on the inlet fitting.

- Place the hose barb fitting on top of the centering ring.
- Secure the fitting with the clamp.


- On rotary vane pumps only, attach the exhaust filter to the exhaust port. Some pumps have filters that screw in. Others require that a centering ring is placed on the fitting and then a clamp secures it similar to the inlet fitting.

- NOTE: Rotary vane vacuum pumps must have an exhaust filter to prevent oil mist from exiting the pump and depleting the oil supply. Fill the pump with oil so the level is between the lines on the sight gauge on the end of the pump housing.
- Verify that the vacuum pump voltage setting agrees with the Concentrator voltage shown on the serial number label.
- Turn the vacuum pump power switch ON. (Rotary Vane or Diaphragm pumps)

- For rotary vane pumps, set the gas ballast to closed or minimum position.

Step 5. Electrical Connections

- Vacuum Pump – The vacuum pump is powered through the CentriVac's controls. If not equipped, attach a reverse IEC connector to the power cord of the vacuum pump so it can plug into the back of the concentrator.
- CentriVac Concentrator – Plug the CentriVac into the electrical outlet that matches the rating shown on the label (lower-rear).

Testing or Start Up

Step 6. Testing the System

- Turn the CentriVac power switch ON.
- Close the lid on the CentriVac
- Close the lid on the Cold Trap (if equipped)
- Press the Run/Stop button
- Once the CentriVac attains speed the vacuum pump should start automatically
- To stop the CentriVac and vent, press the Run/Stop button
- If the CentriVac or Cold Trap fails to perform as expected, consult the User's Manuals for troubleshooting information or call Labconco.

LABCONCO CORPORATION 8811 Prospect Ave, Kansas City, MO 64132
(816) 333-8811 or (800) 821-5525 • (816) 363-0130 fax • labconco.com

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