


Front View



Side View

## # 8032500, 2-Basket Glassware Carts

- Safely transport clean and dirty labware.
- Come with two large or four small interchangeable vinyl-coated wire baskets and removable drip pan below the baskets to catch spills.
- Include 4" diameter casters (two with toe locks). Frame is welded 1" square mar-resistant, powder-coated tubular steel.
- Overall dimensions: 34.9" w x 19.0" d x 36.4" h.
- Fully assembled.
- Lifetime warranty.
- Shipping weight 43 lbs. (20 kg)

Cart  
for:



LABCONCO CORPORATION • KANSAS CITY MO USA 64132  
800-821-5525 / 816-333-8811 • labconco.com

Date:  
09-12-2022

Drawing #  
8032500

Drawn By:  
C. Alumbaugh

SCALE: Scaled to fit