

	FRONT	SIDE	BACK	TOP	
99010-00					<p>Protector Acid Storage Cabinet</p> <ul style="list-style-type: none"> * Supports all Labconco benchtop fume hoods. * Work surface is required. * Durable epoxy-coated steel construction with corrosion resistant polyethylene lined interior. * Attractive glacier white exterior to match Labconco fume hoods. * Manual-closing, non-locking door(s). * Includes vacuum formed PVC liner tray. * Includes epoxy-coated steel interior shelf. * Includes two rear vent connection cutouts for venting cabinet into a fume hood. * Include four leveling feet and one 8" filter panel to increase cabinet depth from 22" to 30". * Work surface, interior shelf and vent kits are sold separately.
99011-00					<p>9901000 - 9901600</p> <ul style="list-style-type: none"> * Standard Height
99012-00					<p>9905000 - 9905300</p> <ul style="list-style-type: none"> * ADA Height
99013-00 (ADA)					
99014-00 (ADA)					

CABINETS
for:



No.	Date

Date:	07-09-2008
Drawing #	9905000 - 9905300
Drawn By:	Luke Savage
Checked By:	
Approved By:	
Sheet No.	