



Front View

Side View

# 8000300, Carboy Caddy

- Safe way to move large carboys or cartons
- Helps eliminate heavy lifting and back strain
- 3" diameter casters provide low profile, elevating shelf only 4.25" above the floor
- Shelf is made of phenolic board with melamine high-pressure laminated surface
- Handle measure 19" w x 38.5" h
- Shelf dimensions: 20" w x 23.6" l
- Supports loads up to 360 lbs.

Cart  
for:



LABCONCO CORPORATION • KANSAS CITY MO USA 64132  
800-821-5525 / 816-333-8811 • labconco.com

Date:  
09-26-2022

Drawing #  
8000300

Drawn By:  
C. Alumbaugh

SCALE:                      Scaled to fit