


Front View



Side View

8047500, Pan Cart

- Two 2" deep powder-coated steel, removable watertight pans transport wet or soiled glassware and other labware.
- Removable rubber mats line pans to prevent loads from shifting
- 4" diameter casters, two with toe locks
- Supports loads up to 400 lbs.
- Dimensions: 34.9" w x 19.0" d x 36.4" h

Cart
for:



LABCONCO CORPORATION • KANSAS CITY MO USA 64132
800-821-5525 / 816-333-8811 • labconco.com

Date:
09-26-2022

Drawing #
8047500

Drawn By:
C. Alumbaugh

SCALE:

Scaled to fit