

Purifier[®] Axiom:[™] Get it All.

LIFETIME SAVINGS
OF MORE THAN

50%*



*Protecting your
laboratory environment*

LABCONCO[®]

*When compared to a Type B2 cabinet

SAVE ENERGY



SAVE AIR



SAVE MONEY



Type A biosafety cabinets recirculate filtered air back into laboratories, and Type B biosafety cabinets are ducted to the outside. The innovative new Axiom, a Type C1 enclosure, is unique in its versatility because it can function in either mode, providing flexibility AND savings.

HERE'S HOW IT WORKS

A special work surface design allows uncontaminated air to be filtered and recirculated when the cabinet is in the Type B mode, saving energy when ducted. Reduced airflow demand equals staggering savings.

THERE'S MORE

The Axiom can be connected to existing building exhaust systems when functioning in Type B mode, eliminating the costs of consulting services, utility space, heavy ducting, exhaust blowers, installation, wiring and roof penetrations. The savings are substantial.

THE BOTTOM LINE

The Axiom is less expensive than a conventional Type B cabinet, more adaptable than a Type A cabinet and can do the work of both.



Watch the video at axiom.labconco.com



Flex Your Options.



Not Your Budget.

- The Axiom Type C1 is the most flexible BSC available. No other BSC can convert from Type A to Type B mode. No need to purchase an entirely new BSC when applications change.
- In Type A mode, the Axiom Type C1 offers the same energy efficiency of a Type A BSC that is not ducted or canopy ducted.
- In Type B mode, the Axiom Type C1 is significantly less expensive to operate than a Type B.

Watch the video at
axiom.labconco.com



*Protecting your
laboratory environment*

LABCONCO®

Labconco Corporation • Kansas City, MO 64132
800.821.5525 or 816.333.8811
www.axiom.labconco.com

2-28-9/12/16