









FRONT



SIDE



BACK



TOP

Protector Solvent Storage Cabinet

- * Supports all Labconco benchtop fume hoods.
- * Work surface is required.
- * Durable dual wall epoxy-coated steel construction with 1.5" air gap.
- * Attractive glider wheel exterior to match Labconco fume hoods.
- * Manual-closing or self-closing door(s) with 3-position locking handle.
- * Meets OSHA 29 CFR 1910, section 1910.106 specifications and NF 97.
- * Includes epoxy-coated steel interior shelf.
- * Includes two rear vent connections with flame arrestors and closure plugs.
- * Includes four leveling feet and one 8" filler panel to increase cabinet depth from 22" to 30".
- * Work surface is sold separately.

9902000 - 9902400

- * Manual-closing door(s) with 3-position locking handle.

9903000 - 9903400

- * Self-closing door(s) with 3-position locking handle

9906000 - 9906100
9906200 - 9906300

- ADA Height

