



## XPert® Balance Enclosures and Bulk Powder Enclosures

### ACCESSORIES



#### XPert Filtered System and Station Replacement HEPA Filters and Accessories\*

Catalog #	Description	For use with:	Shipping Weight
3707900	Replacement HEPA Filter, 99.99% efficient on particles 0.3 micron in size, 18.0" w x 18.0" d x 3.3" high	2-Ft XPert System, Station	8 lbs. (4 kg)
3707901	Replacement HEPA Filter, 99.99% efficient on particles 0.3 micron in size, 30.0" x 18.0" x 3.3" high	3-Ft XPert System, Station 6-Ft XPert System, Station (2 req.)	15 lbs. (7 kg)
3707902	Replacement HEPA Filter, 99.99% efficient on particles 0.3 micron in size, 42.0" w x 18.0" d x 3.3" high	4-Ft XPert System, Station	20 lbs. (9 kg)
3707903	Replacement HEPA Filter, 99.99% efficient on particles 0.3 micron in size 24.0" w x 18.0" d x 3.3" high	5-Ft XPert System, Station (2 required)	13 lbs. (6 kg)
3776003	Replacement Bag for Bag-in/Bag-out, HEPA Filter, 48" circumference x 54" long	2-Ft XPert System, Station	5 lbs. (2 kg)
3776004	Replacement Bag for Bag-in/Bag-out, HEPA Filter, 73" circumference x 96" long	3-Ft XPert System, Station 6-Ft XPert System, Station (2 req.)	7 lbs. (3 kg)
3776005	Replacement Bag for Bag-in/Bag-out, HEPA Filter, 96" circumference x 96" long	4-Ft XPert System, Station	9 lbs. (4 kg)
3776006	Replacement Bag for Bag-in/Bag-out, HEPA Filter, 60" circumference x 72" long	5-Ft XPert System, Station (2 required)	6 lbs. (3 kg)

\* HEPA filter removal and replacement should only be performed by a certified technician.



#### Fluorescent Lights for XPert Enclosures

Rest atop XPert Enclosures for enhanced illumination of the interior. Include auxiliary electrical receptacle with switch for connection of accessories such as a printer or airflow monitor. More than one may be used on larger enclosures. 18.3" w x 10.1" d x 3.0" high.

Catalog #	Description	Shipping Wt. lbs./kg
3909200	115 volts, 60 Hz, 10 amps. Includes IEC, three-wire cord with plug.	15/7
3909201*	230 volts, 50 Hz, 10 amps. Includes IEC, three-wire cord. Plug is required (not included).	15/7



#### 3925000 Utility Shelves

For use with any XPert Enclosure, Station or System. Three epoxy-coated steel shelves attach to the rear baffle perforations to hold a

printer, utensils and other miscellaneous weighing supplies. The utensil shelf has a scalloped front edge to support spatulas, brushes and other items. Shipping weight 5 lbs. (2 kg)



#### 3927700 Bottle Holder

For use with any XPert Enclosure, Station or System. Epoxy-coated steel. Attaches to the rear baffle perforations to hold a squirt bottle or other round container. 3" diameter opening. Shipping weight 1 lb. (0.5 kg)



#### 3927800 Tissue Holder

For use with XPert Enclosure, System or Station. Epoxy-coated steel. Attaches to the rear baffle perforations to hold tissue wipers in the upside down position.

An opening in the holder permits wipers to be dispensed from the bottom. Shipping weight 1 lb. (0.5 kg)



#### 3744000 Ergonomic Chair with Armrests

Chair has 6-way articulating seat and back control for personalized adjustment. Pneumatic mechanism adjusts seat height from 18.25" to 25.75". Five-leg black reinforced composite base rests on 2" ball-bearing casters. Aluminum support ring. Removable arm rests. Black vinyl upholstery. Shipping weight 35 lbs. (16 kg)



#### 3746000 Adjustable Footrest

Elevates feet and permits angle repositioning while in use. 18.5" w x 11.5" d x 8" high. Shipping weight 6 lbs. (3 kg)