



### Replacement ULPA filters and Accessories

Catalog #	Description	For use with:	Shipping Weight
<b>3885700</b>	Replacement ULPA Filter, 99.999% efficient on particles 0.12 micron in size, 18.0" w x 18.0" d x 3.3" high (46 x 46 x 8 cm)	2' XPert Nano Enclosure	8 lbs. (4 kg)
<b>3885701</b>	Replacement ULPA Filter, 99.999% efficient on particles 0.12 micron in size, 30.0" x 18.0" x 3.3" high (76 x 46 x 8 cm)	3' XPert Nano Enclosure 6' XPert Nano Enclosure (2 req.)	15 lbs. (7 kg)
<b>3885702</b>	Replacement ULPA Filter, 99.999% efficient on particles 0.12 micron in size, 42.0" w x 18.0" d x 3.3" high (107 x 46 x 8 cm)	4' XPert Nano Enclosure	20 lbs. (9 kg)
<b>3885703</b>	Replacement ULPA Filter, 99.999% efficient on particles 0.12 micron in size 24.0" w x 18.0" d x 3.3" high (61 x 46 x 8 cm)	5' XPert Nano Enclosure (2 required)	13 lbs. (6 kg)
<b>3776003</b>	Replacement Bag for Bag-in/Bag-out, ULPA Filter, 48" circumference x 54" long (122 x 137 cm)	2' XPert Nano Enclosure	5 lbs. (2 kg)
<b>3776004</b>	Replacement Bag for Bag-in/Bag-out, ULPA Filter, 73" circumference x 96" long (185 x 244 cm)	3' XPert Nano Enclosure 6' XPert Nano Enclosure (2 req.)	7 lbs. (3 kg)
<b>3776005</b>	Replacement Bag for Bag-in/Bag-out, ULPA Filter, 96" circumference x 96" long (244 x 244 cm)	4' XPert Nano Enclosure	9 lbs. (4 kg)
<b>3776006</b>	Replacement Bag for Bag-in/Bag-out, ULPA Filter, 60" circumference x 72" long (152 x 183 cm)	5' XPert Nano Enclosure (2 required)	6 lbs. (3 kg)



### 3744000 Ergonomic Chair with Armrests


Free-floating articulating seat control includes pneumatic seat height adjustment of 19 - 26" (48 - 66 cm) and backrest height adjustment. Seat and backrest tilt independently with separate levers and can free-float or lock into place. Five-leg black composite base with dual wheel casters that resist movement when the chair is unoccupied. Includes height adjustable, self-skinned urethane arm-

rests and 20" (51 cm) diameter adjustable footring. Black vinyl upholstery. Shipping weight 35 lbs. (16 kg)



### 3746000 Adjustable Footrest


Elevates feet and permits angle repositioning while in use. 18.5" w x 11.5" d x 8.0" high (47 x 29 x 20 cm). Shipping weight 6 lbs. (3 kg)



### Canopy Connections

Allow the XPert Nano Enclosure to be thimble-ducted to the outside. The canopy's air gap provides a means of balancing air by allowing room air to enter the ductwork and mix with the filtered air

being exhausted from the enclosure. **House exhaust or remote blower and ductwork are required (not included). Requires installation.**



### 5234600 Balance Vibration Isolator

Minimizes vibration to sensitive balances. Includes a 13.4" w x 20.4" d x 2.7" thick (34 x 52 x 7 cm) marble slab, four isolator pads, and a type 304 stainless steel protective cover. Shipping weight 75 lbs. (34 kg)

Catalog #	For use with:	Collar Dia.	Shipping Weight
<b>3924400</b>	2' XPert	6" (15 cm)	15 lbs. (7 kg)
<b>3924401</b>	3' XPert	6" (15 cm)	18 lbs. (8 kg)
<b>3924402</b>	4' XPert	6" (15 cm)	21 lbs. (10 kg)
<b>3924403</b>	5' XPert	10" (25 cm)	25 lbs. (11 kg)
<b>3924404</b>	6' XPert	10" (25 cm)	30 lbs. (14 kg)