


Zero Pressure Weathercap,
12” Diameter
Model: # 5622100

Description: PVC cap adds little static pressure to the exhaust system. Allows for vertical discharge of effluent air for dispersion away from the building.

Equivalent Resistance of
Straight Duct: 5



Zero Pressure
Weathercap for:



LABCONCO CORPORATION • KANSAS CITY MO USA 64132
800-821-5525 / 816-333-8811 • labconco.com

Date:
07-12-2024

Drawing #
5622100

Drawn By:
C. Alumbaugh